INSTRUCTIONS FOR USE OF AB DENTAL GERMANY UG REUSEABLE INSTRUMENTS

DISCLAIMER: AB Dental Germany UG products are intended for use only by certified dentists and authorized personnel with specific implant training.

GENERAL

These instruction for use for AB Dental Germany UG instruments products includes the following categories: Tools, Drills and Kit Boxes. In this IFU the term instruments will describe all these products. Instruments are reusable and should be cleaned, disinfected, and sterilized prior to any use.

INDICATIONS

AB Dental Germany UG instruments are intended for use in surgery and restorative applications in order to restore the patient's chewing function. Specific indication per each instrument category is described below.

Patient Target

The products are intended for women or men over the age of 18 years old.

WARNINGS

Use of non-sterile components may lead to infection of tissues infectious diseases.

CAUTION

- Federal (USA) law restricts the sale of this device to, or on the order of a licensed physician or dentist.
- All the instruments that are used during the clinical procedure or technician laboratory process, must be maintained in good condition and the user should take care to their physical condition to avoid any damage to other associated products such as implants and prosthetic parts or damage to the bone.
- Worn instruments must be disposed/discarded according to local authority regulations and environmental requirements or returned to AB Dental Germany UG.
- Always wear gloves when handling contaminated instruments.
- Eye protection must be worn to protect against ejected particles.
- Surgical mask must be worn to avoid inhalation of aerosol or generated dust.
- Ensure that the drills or relevant tools are fully seated and gripped in the hand piece collet before use.
- Maintain hand piece in good working order and according to manufacturer instructions.
- Do not exceed the maximum drill speed shown in the recommended drilling protocol in products catalog, to prevent overheating and associated complications.
- Avoid excessive drilling duration to prevent overheating and associated complications.
- Move the drill continuously when in use to avoid localized heating.
- We advise to replace the drill every 10 uses.
- Avoid cleaning, disinfecting and sterilizing together, tools and drills from different materials.
- Avoid contact with solutions that main ingredient has high concentration of: phenol alcohol, chlorine, aldehyde, acid, quaternary ammonia, hydrogen peroxide, saline solution and water with high concentration of chlorine.
- Use only specifically formulated cleaning and disinfectant agents.

CONTRAINDICATIONS

Customary observations should be made of the contraindications associated with tools and drills materials used in oral surgery. The patient's general health and suitability for oral surgery must be assessed by the general practitioner.

STORAGE AND HANDLING

The products must be inspected prior to usage. The products should be stored in clean environment, away from direct sunlight. Incorrect storage might damage the products characteristics and lead to product failure.

CLEANING AND STERILIZATION PROCESS

CLEANING AND DISINFECTION INSTRUCTIONS:

1. Directly after use (within 30 min.), disassemble multi-piece instruments, wipe and rinse the instrument to remove gross soil, blood and tissue before cleaning.
Use a nylon brush to remove debris from instruments, particularly if caught in blind hole or lumen and flush with a hypodermic needle.
2. Immerse the drills and/or tools in a mild, pH-neutral commercially available detergent for use with dental instrument according to manufacturer's instructions.
Use a nylon brush to remove debris from instruments, particularly if caught in blind hole or lumen and flush with a hypodermic needle.
Rinse the instruments under a hard stream of water for a minimum duration of 3 minutes.
3. Immerse the drills and/or tools in a mild, pH-neutral commercially available disinfectant for use with dental instrument according to manufacturer's instruction and flush with a hypodermic needle.
Rinse the instruments under a hard stream of water for a minimum duration of 3 minutes.
4. Dry the instruments with clean lint-free single use wipes or dry heat not exceeding 134°C/275°F.

STERILIZATION INSTRUCTIONS:

5. Place the instruments for sterilization inside commercially available autoclave pouch.

6. Sterilize the instruments by steam per the autoclave's instructions. Set the Autoclave to 134C/275°F for 6 minutes followed by a 30 minutes' dry cycle. Distilled water should be used to avoid surface stains. Make sure before use that the elements, inside the autoclave, are not rusted.
7. Chemclave is NOT recommended.

PRODUCT FAMILIES – INTENDED USE AND PRODUCT DESCRIPTION

Drills

The drills are intended to prepare the osteotomy in the jawbone for insertion of the dental implant. The drills are made from medical grade stainless-steel with or without coating. The drills are connected to handpiece. The drills have been marked with sizes and color coded for ease of identification. Size marking and color coding should be used to select proper drill for each procedure. Follow the recommended drill protocol according to the implant characteristics.

Following are all families of drills:

Product family	Description	Intended Use
TMD	Marking drills	Mark the drilling site
TPD TD	Pilot drill, Straight drills	Drill to desired/implant length according to depth marks
TPDD TDD	Pilot coated drill Straight coated drills	Drill to desired/implant length according to depth marks
TDS	Step coated drills	Drill to desired/implant length according to depth marks
TDCSt	Conical stopper drills	Tapered/Conical shape drill with dedicated length according to implant length
TDG	Guide surgery drill	Drill through guide stent according to desired/implant length
TDCS	Countersink drills	Drill cortical bone
TDTI	Trephine drills	Drill around implant to remove the implant

Tools

Following are all families of tools:

Product family or specific catalog number/numbers	Description	Intended Use
T3, T3G	Ratchet implant drivers	Thread the implant to the osteotomy
T5	Contra angle (rotary motor mount) implant or prosthetic drivers	Thread the implant to the osteotomy or thread the prosthetic products to implant or other products
T1, T2, T13	Prosthetic drivers	Thread the prosthetic products to implant or other products
T8, T8c	Ratchet and torque ratchet	Transfer moment torque to ratchet tools
T9	Depth guide	Probe to view the depth of osteotomy
T10	Handle	Transfer moment torque to ratchet tools
T11	Mallet	To strike on other tools that will help to break bone
T15	Broken Screw Removal Tool	Extract broken screw
T16	Implant angle position	Visual aid for angulation
T17	Tissue punch	Cut Gingival tissue
T18	Implant Removal Tool	Extract implant with worn Hex
T22	Abutment gripper	Technician too to hold an abutment and customize it
T4	Retrieving screws	Extract stuck abutment form implant
TP	Pins	Hold guide in place
TH	Handles	Guide adapter for sleeves and drills
TDE	Drill extension	Contra angle (rotary) adapter/extender
D2, D20, D3, PK-D2	Impression coping Transfers	Close and open tray impression taking aids (see additional details below)
D1	Analog/Digital Model Analog	Implant and prosthetics parts replicas for models
D4	Prosthetic Kit	Sleeve for Multiunit (P64, P14C-4 or P16C-4) + Close tray impression taking aids (see additional details below)
P3,SC/P14,SC/P64,SC	Scan Abutments/Bodies	Transfer the implant or screw retained abutment position and orientation to a digital file for later prosthetic restoration planning (see additional details below)

Prosthetic Kits D4 and PK

D4-3.75,X; D4-3,X; PK-D2

The D4 prosthetic kits include dedicated abutment for implants, three plastic transfer PK-D2 (a disposable plastic snap transfer) and implant analog.

D4-P64; D4-P14C; PK-D2

The D4 and PK prosthetic kits include dedicated sleeve (not for final restoration) for screw retained abutments and three plastic transfer PK-D2.

Following the connection of the dedicated abutment or dedicated sleeve the PK-D2 is snapped on the relevant product and allows the user to take close tray impression. The relevant product/analog is snap to the PK-D2 plastic transfers, that are in the impression and then sent to the technician lab allowing the lab to create cast model.

It is recommended to sterilize the metal parts prior to the use.

The plastic transfers PK-D2 and PK-P2 should not be sterilized.

Scan Abutments/Bodies

P3,SC/P14,SC/P64,SC

The scan abutments/bodies intended to transfer the implant/Multiunit (screw retained abutment) position and orientation data (in patient mouth or cast model) to a digital file for later prosthetic restoration planning.

The scan abutments/bodies can be used as intraoral scan abutment or cast model scan abutment/bodies and intended for multiple use.

The scan abutments/bodies are made of Biocompatible PEEK with at least 2 flat surfaces for easy and accurate alignment/registration of scan image with the library model. The scan abutment/body correlate with the implant or screw retained connection. For better scanning precision we recommend locating the flat surfaces in palatal/lingual orientation. For use with angled implant/abutment, the short flat face on the scan abutment should be aligned with the slanted (Buccal) side of the base.

Scan abutment/body is secured with screw to the implant or screw retained abutment with the recommended closure torque according to the relevant platform.

The different scan abutments/bodies correspond with the respective platform (implants or screw retained) and differ in shapes/sizes. The scanning process should be performed as recommended by scanning/CADCAM system. It is important to choose the correct platform in software and corresponding type for the chosen restoration (engaging/non-engaging). After use the scan abutment can be loosened and placed gently on the storage.

Scan abutments must be cleaned and sterile prior to use in patient mouth. Please follow cleaning/sterile instructions.

Kit boxes

The Kit boxes are intended to organize and contain the dental instruments during sterilization process, surgery, storage and transport.

DISPOSAL

No longer usable parts must be disposed/discarded according to local authority regulations and environmental requirements. In case contaminated products need to be returned to AB Dental Germany UG, please follow the company return policy.

NOTE

Issues that arise in relation to the products should be reported together with the related product to the AB Dental Germany UG representative.

In case that any serious incident occurs in relation with the products, the user and/or patient should report to AB Dental Germany UG and to the competent authority of the Member State in which the serious incident occurs.

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CE

CE Mark for
Class I Devices

	Caution Indicated the need to consult the instructions for use for important cautionary
	CE marking of conformity
	Batch code Lot number
	Catalogue number
	Consult instructions for use: www.ab-dent.com/ifu
	Date of manufacture
	Manufacturer
	Non-sterile
	Subject to prescription
	Medical Device
	Keep away from sunlight
	Do not use if package is damaged

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Ürünlerle ilgili olarak ortaya çıkan sorunlar, AB Dental Germany UG temsilcisine bildirilmelidir.
AB Dental Germany UG (Haftungsbeschränkt)