






Superior Implant Technology

I55

A NEW IMPLANT WITH
THE OPTIMAL DESIGN FOR
IMMEDIATE LOADING





I55

A NEW IMPLANT WITH THE OPTIMAL DESIGN FOR IMMEDIATE LOADING

Features & Benefits:

- ▀ Two separate external threads:
 - Flat thread that enables the strongest initial stability
 - Sharp thread that enables the insertion of the implant easily and with minimal trauma to the bone
- ▀ Narrow apex allowing implant insertion after drilling only with the pilot drill
- ▀ Large surface area facilitating the osseointegration process

Sizes:

Cat no.	Diameter	Length
I55	3	10, 11.5, 13, 16
I55	3.5	10, 11.5, 13, 16
I55	3.75	8, 10, 11.5, 13, 16
I55	4.2	8, 10, 11.5, 13, 16
I55	4.5	6, 8, 10, 11.5, 13, 16
I55	5	6, 8, 10, 11.5, 13, 16

www.ab-dent.com