



To be the leading innovative authority protecting and advancing national health and safety

PUBLIC ENQUIRY - SINGAPORE MEDICAL DEVICE REGISTER (SMDR)

Medical Device Device Category Registrant Product Owner Importer Wholesaler Advanced Search

Device Info

Device Name: [A.B. Dental Devices A.B. Dental Implant System](#) [AB DENTAL DEVICES LTD]

Description: A.B. Dental Devices' implants are designed for prosthodontic reconstruction. They are to be inserted into areas of the maxillary or mandibular arch locations, for completely or partially edentulous patients. Patients must have adequate alveolar bone width to support the width and length of the implant (minimum 1 mm circumferential and 2 mm apical). Sufficient volume and quality of bone used for placement of implant may be determined through adequate radiological data used by implant dentistry. For use in the following cases: single tooth loss, free-end saddle replacement, total tooth loss and stabilization of removable appliances. A thorough evaluation of the following must be made with their relation to the implant site: vital blood vessel location, nerve location, maxillary sinus and soft tissue spaces. Implants in the maxilla should not puncture the sinus, unless sinus augmentation is performed.

Medical Specialty Area: Dental

HS Code: 90212900

HSA Product code: HSAMDC00103

Medical Device Class: Class C medical device

Device Registration No: DE0016409

Registration Date: 13/10/2014

Change Notification Approval Date: Not Applicable

Expiry Date: 12/10/2015

Device System Info: [Sterile Medical Device], Dental implants are sterilized by gamma radiation.

Product Owner

1. [AB DENTAL DEVICES LTD](#) [AB DENTAL DEVICES LTD] 3 HABOSEM STREET, MINRAV BUILDINGS, NORTH INDUSTRIAL AREA, ASHDOD 77...

Registrant

1. [CHEE SANG DENTAL SUPPLY PTE LTD](#) 57 UBI AVENUE 1, UBI CENTRE, #04-01, SINGAPORE 408936

Imported By

1. [CHEE SANG DENTAL SUPPLY PTE LTD](#) [ES0000892] 57 UBI AVENUE 1, UBI CENTRE, #04-01, SINGAPORE 408936

Models

Product Code for Accessories:

No.	Model Name	Identifier
1	Screw Type Implant	I2-3.25,8
2	Screw Type Implant	I2-3.25,10
3	Screw Type Implant	I2-3.25,11.5
4	Screw Type Implant	I2-3.25,13
5	Screw Type Implant	I2-3.25,16
6	Screw Type Implant	I2-3.75,8
7	Screw Type Implant	I2-3.75,10
8	Screw Type Implant	I2-3.75,11.5

9	Screw Type Implant	I2-3.75,13
10	Screw Type Implant	I2-3.75,16
11	Screw Type Implant	I2-4.2,8
12	Screw Type Implant	I2-4.2,10
13	Screw Type Implant	I2-4.2,11.5
14	Screw Type Implant	I2-4.2,13
15	Screw Type Implant	I2-4.2,16
16	Screw Type Implant	I2-4.2,18
17	Screw Type Implant	I2-4.2,20
18	Screw Type Implant	I2-4.5,8
19	Screw Type Implant	I2-4.5,10
20	Screw Type Implant	I2-4.5,11.5
21	Screw Type Implant	I2-4.5,13
22	Screw Type Implant	I2-4.5,16
23	Screw Type Implant	I2-5,8
24	Screw Type Implant	I2-5,10
25	Screw Type Implant	I2-5,11.5
26	Screw Type Implant	I2-5,13
27	Screw Type Implant	I2-5,16
28	Screw Type Implant	I2-6,8
29	Screw Type Implant	I2-6,10
30	Screw Type Implant	I2-6,11.5
31	Screw Type Implant	I2-6,13
32	Screw Type Implant	I2-6,16
33	Conical Implant	I5-3,10
34	Conical Implant	I5-3,11.5
35	Conical Implant	I5-3,13
36	Conical Implant	I5-3,16
37	Conical Implant	I5-3.2,8
38	Conical Implant	I5-3.2,10
39	Conical Implant	I5-3.2,11.5
40	Conical Implant	I5-3.2,13
41	Conical Implant	I5-3.2,16
42	Conical Implant	I5-3.75,8
43	Conical Implant	I5-3.75,10
44	Conical Implant	I5-3.75,11.5
45	Conical Implant	I5-3.75,13
46	Conical Implant	I5-3.75,16
47	Conical Implant	I5-4.2,8
48	Conical Implant	I5-4.2,10
49	Conical Implant	I5-4.2,11.5
50	Conical Implant	I5-4.2,13
51	Conical Implant	I5-4.2,16
52	Conical Implant	I5-4.5,6
53	Conical Implant	I5-4.5,8
54	Conical Implant	I5-4.5,10
55	Conical Implant	I5-4.5,11.5
56	Conical Implant	I5-4.5,13
57	Conical Implant	I5-4.5,16
58	Conical Implant	I5-5,6
59	Conical Implant	I5-5,8
60	Conical Implant	I5-5,10

61	Conical Implant	I5-5,11.5
62	Conical Implant	I5-5,13
63	Conical Implant	I5-5,16
64	Conical Implant	I5-6,6
65	Conical Implant	I5-6,8
66	Conical Implant	I5-6,10
67	Conical Implant	I5-6,11.5
68	Conical Implant	I5-6,13
69	Conical Implant	I5-6,16
70	Conical Implant	I5-7,6
71	Conical Implant	I5-7,8
72	Conical Implant	I5-7,10
73	Conical Implant	I5-7,11.5
74	Conical Implant	I5-7,13
75	Conical Implant	I5-8,6
76	Conical Implant	I5-8,8
77	Conical Implant	I5-8,10
78	Conical Implant	I5-8,11.5
79	Conical Implant	I5-8,13
80	Conical Implant	I5-9,6
81	Conical Implant	I5-9,8
82	Conical Implant	I5-9,10
83	Narrow Integral Implant	I6-2.4,10
84	Narrow Integral Implant	I6-2.4,11.5
85	Narrow Integral Implant	I6-2.4,13
86	Narrow Integral Implant	I6-2.4,16
87	Narrow Integral Implant	I6-3,8
88	Narrow Integral Implant	I6-3,10
89	Narrow Integral Implant	I6-3,11.5
90	Narrow Integral Implant	I6-3,13
91	Narrow Integral Implant	I6-3,16
92	Narrow Integral Implant	I6-3.2,8
93	Narrow Integral Implant	I6-3.2,10
94	Narrow Integral Implant	I6-3.2,11.5
95	Narrow Integral Implant	I6-3.2,13
96	Narrow Integral Implant	I6-3.2,16
97	One Piece Ball Attachment Implant	I6b-2.4,10
98	One Piece Ball Attachment Implant	I6b-2.4,11.5
99	One Piece Ball Attachment Implant	I6b-2.4,13
100	One Piece Ball Attachment Implant	I6b-2.4,16
101	Narrow Implant	I6B-3,8
102	Narrow Implant	I6B-3,10
103	Narrow Implant	I6B-3,11.5
104	Narrow Implant	I6B-3,13
105	Narrow Implant	I6B-3,16
106	Narrow Implant	I6B-3.2,8
107	Narrow Implant	I6B-3.2,10
108	Narrow Implant	I6B-3.2,11.5
109	Narrow Implant	I6B-3.2,13
110	Narrow Implant	I6B-3.2,16
111	Narrow Implant	I6BI-3.2,8
112	Narrow Implant	I6BI-3.2,10

113	Narrow Implant	I6BI-3.2,11.5
114	Narrow Implant	I6BI-3.2,13
115	Narrow Implant	I6BI-3.2,16
116	Narrow Implant	I6BI-3,8
117	Narrow Implant	I6BI-3,10
118	Narrow Implant	I6BI-3,11.5
119	Narrow Implant	I6BI-3,13
120	Narrow Implant	I6BI-3,16
121	Integral Implant	I7-3,8
122	Integral Implant	I7-3,10
123	Integral Implant	I7-3,11.5
124	Integral Implant	I7-3,13
125	Integral Implant	I7-3,16
126	Integral Implant	I7-3.75,6
127	Integral Implant	I7-3.75,8
128	Integral Implant	I7-3.75,10
129	Integral Implant	I7-3.75,11.5
130	Integral Implant	I7-3.75,13
131	Integral Implant	I7-3.75,16
132	Integral Implant	I7-4.2,6
133	Integral Implant	I7-4.2,8
134	Integral Implant	I7-4.2,10
135	Integral Implant	I7-4.2,11.5
136	Integral Implant	I7-4.2,13
137	Integral Implant	I7-4.2,16
138	Integral Implant	I7-5,6
139	Integral Implant	I7-5,8
140	Integral Implant	I7-5,10
141	Integral Implant	I7-5,11.5
142	Integral Implant	I7-5,13
143	Integral Implant	I7-5,16
144	Integral Implant	I7-6,6
145	Integral Implant	I7-6,8
146	Integral Implant	I7-6,10
147	Integral Implant	I7-6,11.5
148	Integral Implant	I7-6,13
149	Integral Implant	I7-6,16
150	Trapaze Implant	I10-3.75,8
151	Trapaze Implant	I10-3.75,10
152	Trapaze Implant	I10-3.75,11.5
153	Trapaze Implant	I10-3.75,13
154	Trapaze Implant	I10-3.75,16
155	Trapaze Implant	I10-4.2,8
156	Trapaze Implant	I10-4.2,10
157	Trapaze Implant	I10-4.2,11.5
158	Trapaze Implant	I10-4.2,13
159	Trapaze Implant	I10-4.2,16
160	Trapaze Implant	I10-4.5,8
161	Trapaze Implant	I10-4.5,10
162	Trapaze Implant	I10-4.5,11,5
163	Trapaze Implant	I10-4.5,13
164	Trapaze Implant	I10-4.5,16

165	Trapeze Implant	I10-5,8
166	Trapeze Implant	I10-5,10
167	Trapeze Implant	I10-5,11.5
168	Trapeze Implant	I10-5,13
169	Trapeze Implant	I10-5,16
170	Short and Wide Implant	I15-5,5
171	Short and Wide Implant	I15-6,4
172	Short and Wide Implant	I15-7,4
173	Short and Wide Implant	I15-8,4
174	Short and Wide Implant	I15-9,4
175	Short and Wide Implant	I15-9,5
176	Short and Wide Implant	I15-6,5
177	Short and Wide Implant	I15-7,5
178	Short and Wide Implant	I15-8,5
179	Healing Cap for I6-B Implant	P0-I6B
180	Healing Cap for Aesthetic Abutment with Hex	P0-P6H
181	Titanium Healing Cap	P0-3,2
182	Titanium Healing Cap	P0-3,3
183	Titanium Healing Cap	P0-3,4
184	Titanium Healing Cap	P0-3,5
185	Titanium Healing Cap	P0-3,6
186	Titanium Healing Cap	P0-3.75,2
187	Titanium Healing Cap	P0-3.75,3
188	Titanium Healing Cap	P0-3.75,4
189	Titanium Healing Cap	P0-3.75,5
190	Titanium Healing Cap	P0-3.75,7
191	Titanium Healing Cap	P0-4.5,3
192	Titanium Healing Cap	P0-4.5,5
193	Titanium Healing Cap	P0-4.5,7
194	Titanium Healing Cap	P0/2-3.75,2
195	Titanium Healing Cap	P0/2-3.75,3
196	Titanium Healing Cap	P0/2-3.75,4
197	Titanium Healing Cap	P0/2-3.75,5
198	Titanium Healing Cap	P0/2-3.75,6
199	Titanium Healing Cap	P0/2-3.75,7
200	Wide Titanium Healing Cap	POW-3.75,3
201	Wide Titanium Healing Cap	POW-3.75,4
202	Wide Titanium Healing Cap	POW-3.75,5
203	Wide Titanium Healing Cap	POW-3.75,6
204	Wide Titanium Healing Cap	POW-4.5,5
205	Wide Titanium Healing Cap	POW/2-3.75,3
206	Wide Titanium Healing Cap	POW/2-3.75,4
207	Wide Titanium Healing Cap	POW/2-3.75,5
208	Wide Titanium Healing Cap	POW/2-3.75,6
209	Narrow Titanium Healing Cap	PON-3.75,2
210	Narrow Titanium Healing Cap	PON-3.75,3
211	Narrow Titanium Healing Cap	PON-3.75,4
212	Narrow Titanium Healing Cap	PON-3.75,5
213	Narrow Titanium Healing Cap	PON-3.75,6
214	Narrow Titanium Healing Cap	PON-3.75,7
215	Healing Cap for-I6B implant	P0-I6B,3
216	Healing Cap for-P14	P0-P14,3

217	Healing Cap for-P14	P0-P14,5
218	Non- Hexed Abutment	P1-3.75,9
219	Non- Hexed Abutment	P1-3.75,13
220	Non- Hexed Abutment	P1-3.75,15
221	Non- Hexed Abutment	P1-4.5,9
222	Non- Hexed Abutment	P1H-3.75,9
223	Non- Hexed Abutment	P1H-3.75,13
224	Non- Hexed Abutment	P1H-3.75,15
225	Plastic Sleeve for I6B Implant	P2-I6B
226	Plastic Sleeve w/Fixation Screw	P2H-4.5,9
227	Plastic Sleeve for I6B	P2H-I6B
228	Plastic Sleeve w/Fixation Screw – Straight w/ Hex	P2H- 3.75,9
229	Plastic Sleeve w/Fixation Screw – Straight	P2- 3.75
230	Plastic Sleeve narrow w/Fixation Screw – Straight	P2N- 3,15
231	Plastic Sleeve w/Fixation Screw – Straight w/ Hex	P2NH-3,15
232	Plastic Sleeve w/Fixation Screw – Straight	P2-4.5,9
233	Plastic Sleeve w/Fixation Screw – Straight w/ Hex	P2NH-3.75,15
234	Plastic Sleeve w/Fixation Screw – Wide w/ Hex	P2WH-3.75,9
235	Plastic Sleeve w/Fixation Screw - Rounded	P2N-3.75,15
236	Straight Anatomic Plastic Sleeve	P2-P3S-3,1
237	Straight Anatomic Plastic Sleeve	P2-P3S-3,2
238	Straight Anatomic Plastic Sleeve	P2-P3S-3,3
239	Angular Anatomic Plastic Sleeve	P2-P4S-3,15,1
240	Angular Anatomic Plastic Sleeve	P2-P4S-3,15,2
241	Angular Anatomic Plastic Sleeve	P2-P4S-3,15,3
242	Angular Anatomic Plastic Sleeve	P2-P4S-3,25,1
243	Angular Anatomic Plastic Sleeve	P2-P4S-3,25,2
244	Angular Anatomic Plastic Sleeve	P2-P4S-3,25,3
245	Straight Anatomic Plastic Sleeve	P2-P3S-3.75,1
246	Straight Anatomic Plastic Sleeve	P2-P3S-3.75,2
247	Straight Anatomic Plastic Sleeve	P2-P3S-3.75,3
248	Angular Anatomic Plastic Sleeve	P2-P4S-3.75,15,1
249	Angular Anatomic Plastic Sleeve	P2-P4S-3.75,15,2
250	Angular Anatomic Plastic Sleeve	P2-P4S-3.75,15,3
251	Angular Anatomic Plastic Sleeve	P2-P4S-3.75,25,1
252	Angular Anatomic Plastic Sleeve	P2-P4S-3.75,25,2
253	Angular Anatomic Plastic Sleeve	P2-P4S-3.75,25,3
254	Plastic Conical Sleeve (White)	P2/HK
255	Anti-rotation Abutment for I6B Implant	P3-I6B
256	Screw for anti-rotation Abutment for I6B Implant	P3-I6Ba
257	Anti-rotation Abutment	P3-3,9
258	Anti-rotation Abutment	P3-3.75,5
259	Anti-rotation Abutment - Perfect Fit	P3-3.75,5-P
260	Anti-rotation Abutment	P3-3.75,7
261	Anti-rotation Abutment - Perfect Fit	P3-3.75,7-P
262	Anti-rotation Abutment	P3-3.75,9
263	Anti-rotation Abutment - Perfect Fit	P3-3.75,9-P
264	Anti-rotation Abutment	P3-3.75,12
265	Anti-rotation Abutment - Perfect Fit	P3-3.75,12-P
266	Anti-rotation Abutment	P3-3.75,15
267	Anti-rotation Abutment with conical platform	P3C-3.75,15
268	Anti-rotation Abutment - Perfect Fit	P3-3.75,15-P

269	Anti-rotation Abutment	P3-4.5,9
270	Anti-rotation Abutment - Perfect Fit	P3-4.5,9-P
271	Anti-rotation Abutment	P3-4.5,12
272	Anti-rotation Abutment - Perfect Fit	P3-4.5,12-P
273	Anti-rotation Abutment	P3-4.5,15
274	Anti-rotation Abutment - Perfect Fit	P3-4.5,15-P
275	Anti-rotation Abutment	P3-5,5
276	Anti-rotation Abutment - Perfect Fit	P3-5,5-P
277	Anti-rotation Abutment	P3-5,7
278	Anti-rotation Abutment - Perfect Fit	P3-5,7-P
279	Anti-rotation Abutment	P3-5,9
280	Anti-rotation Abutment - Perfect Fit	P3-5,9-P
281	Anti-rotation Abutment	P3-6,9
282	Anti-rotation Abutment - Perfect Fit	P3-6,9-P
283	Anatomic Anti-rotation Abutment	P3S-3.75,1
284	Anatomic Anti-rotation Abutment - Perfect Fit	P3S-3.75,1-P
285	Anatomic Anti-rotation Abutment	P3S-3.75,2
286	Anatomic Anti-rotation Abutment - Perfect Fit	P3S-3.75,2-P
287	Anatomic Anti-rotation Abutment	P3S-3.75,3
288	Anatomic Anti-rotation Abutment - Perfect Fit	P3S-3.75,3-P
289	Anatomic Anti-rotation Abutment	P3S-4.5,1
290	Anatomic Anti-rotation Abutment - Perfect Fit	P3S-4.5,1-P
291	Anatomic Anti-rotation Abutment	P3S-4.5,3
292	Anatomic Anti-rotation Abutment - Perfect Fit	P3S-4.5,3-P
293	Temporary PEEK Anatomic Anti-rotation Abutment	P3S-PEEK,1
294	Temporary PEEK Anatomic Anti-rotation Abutment	P3S-PEEK,2
295	Temporary PEEK Anatomic Anti-rotation Abutment	P3S-PEEK,3
296	Anatomic Anti-rotation Abutment	P3K-3.75,1
297	Anatomic Anti-rotation Abutment	P3K-3.75,2
298	Anatomic Anti-rotation Abutment	P3K-3.75,3
299	Anatomic Anti-rotation Abutment	P3K-3.75,4
300	Narrow Anti-rotation Abutment	P3N-3.75,5
301	Narrow Anti-rotation Abutment- Perfect Fit	P3N-3.75,5-P
302	Narrow Anti-rotation Abutment	P3N-3.75,7
303	Narrow Anti-rotation Abutment- Perfect Fit	P3N-3.75,7-P
304	Narrow Anti-rotation Abutment	P3N-3.75,9
305	Narrow Anti-rotation Abutment - Perfect Fit	P3N-3.75,9-P
306	Wide Anti-rotation Abutment	P3W-3.75,9
307	Wide Anti-rotation Abutment - Perfect Fit	P3W-3.75,9-P
308	Wide Anti-rotation Abutment	P3W-3.75,12
309	Wide Anti-rotation Abutment - Perfect Fit	P3W-3.75,12-P
310	Wide Anti-rotation Abutment	P3W-3.75,15
311	Wide Anti-rotation Abutment - Perfect Fit	P3W-3.75,15-P
312	Wide Anti-rotation Abutment	P3W-4.5,9
313	Wide Anti-rotation Abutment - Perfect Fit	P3W-4.5,9-P
314	Wide Anti-rotation Abutment	P3SW-3.75,1
315	Wide Anti-rotation Abutment	P3SW-3.75,2
316	Wide Anti-rotation Abutment	P3SW-3.75,3
317	Zirconium Anatomic Abutment	P3SZ-3,1
318	Zirconium Anatomic Abutment	P3SZ-3,2
319	Zirconium Anatomic Abutment	P3SZ-3,3
320	Zirconium Anatomic Abutment	P3SZ-3.75,1

321	Zirconium Anatomic Abutment	P3SZ-3.75,2
322	Zirconium Anatomic Abutment	P3SZ-3.75,3
323	Abutment screw - Long	P3a-L
324	Abutment screw - Short	P3a-S
325	Angular Abutment 15°	P4-3,15
326	Angular Abutment 25°	P4-3,25
327	Angular Abutment 15°	P4-3.75,15
328	Angular Abutment 15° - Perfect Fit	P4-3.75,15-P
329	Angular Abutment 25° - Perfect Fit	P4-3.75,25-P
330	Angular Abutment 35°	P4-3.75,35
331	Angular Abutment 45°	P4-3.75,45
332	Angular Abutment 15°	P4-4.5,15
333	Angular Abutment 15°	P4-4.5,15-P
334	Angular Abutment 15°	P4-5,15
335	Angular Abutment 15°	P4-6,15
336	Long Angular Abutment 15°	P4L-3.75,15
337	Long Angular Abutment 15° - Perfect Fit	P4L-3.75,15-P
338	Long Angular Abutment 25°	P4L-3.75,25
339	Long Angular Abutment 25° - Perfect Fit	P4L-3.75,25-P
340	Narrow Angular Abutment 15°	P4N-3.75,15
341	Anatomic Angular Abutment for I6B Implant	P4S-I6B
342	Temporary PEEK Anatomic	P4S-Peek 15,1
343	Temporary PEEK Anatomic	P4S-Peek 15,2
344	Temporary PEEK Anatomic	P4S-Peek 15,3
345	Temporary PEEK Anatomic	P4S-Peek 25,1
346	Temporary PEEK Anatomic	P4S-Peek 25,2
347	Temporary PEEK Anatomic	P4S-Peek 25,3
348	Anatomic Angular Abutment 15° with Shoulder	P4S-3.75,15-1
349	Anatomic Angular Abutment 15° with Shoulder-Perfect	P4S-3.75,15-1-P
350	Anatomic Angular Abutment 15° with Shoulder	P4S-3.75,15-2
351	Anatomic Angular Abutment 15° with Shoulder-Perfect	P4S-3.75,15-2-P
352	Anatomic Angular Abutment 15° with Shoulder	P4S-3.75,15-3
353	Anatomic Angular Abutment 15° with Shoulder - Perfect Fit	P4S-3.75,15-3-P
354	Anatomic Angular Abutment 25° with Shoulder	P4S-3.75,25-1
355	Anatomic Angular Abutment 25° with Shoulder	P4S-3.75,25-2
356	Anatomic Angular Abutment 25° with Shoulder	P4S-3.75,25-3
357	Zirconium Anatomic Angular Abutment	P4SZ-3,1
358	Zirconium Anatomic Angular Abutment	P4SZ-3,2
359	Zirconium Anatomic Angular Abutment	P4SZ-3,3
360	Zirconium Anatomic Angular Abutment	P4SZ-3.75,1
361	Zirconium Anatomic Angular Abutment	P4SZ-3.75,2
362	Zirconium Anatomic Angular Abutment	P4SZ-3.75,3
363	Zirconium Anatomic Angular Abutment	P4SZ-3.75,25-1
364	Zirconium Anatomic Angular Abutment	P4SZ-3.75,25-2
365	Zirconium Anatomic Angular Abutment	P4SZ-3.75,25-3
366	Screw for Abutment	P4a
367	Metal Cap	P5a
368	Silicon Cap - soft	P5b-1
369	Silicon Cap- diameter 1.8mm	P5b-1.8
370	Silicon Cap - medium	P5b-2
371	Silicon Cap- diameter 2.2mm	P5b-2.2
372	Silicon Cap - hard	P5b-3

373	Silicon Cap	P5b-t
374	Protective disk	P5d
375	Protective disk	P5pd
376	Ball Attachment Abutment - a+ b+c	P5-3,1set
377	Ball Attachment Abutment - a+ b+c	P5-3,2set
378	Ball Attachment Abutment - a+ b+c	P5-3,3set
379	Ball Attachment Abutment - a+ b+c	P5-3,4set
380	Ball Attachment Abutment - a+ b+c	P5-3,5set
381	Ball Attachment Abutment - a+ b+c	P5-3,6set
382	Ball Attachment Abutment - a+ b+c	P5-3.75,1set
383	Ball Attachment Abutment - a+ b+c	P5-3.75,2set
384	Ball Attachment Abutment - a+ b+c	P5-3.75,3set
385	Ball Attachment Abutment - a+ b+c	P5-3.75,4set
386	Ball Attachment Abutment - a+ b+c	P5-3.75,5set
387	Ball Attachment Abutment - a+b+c	P5-3.75,6set
388	Ball Attachment Abutment	P5-3.75,1
389	Ball Attachment Abutment	P5-3.75,2
390	Ball Attachment Abutment	P5-3.75,3
391	Ball Attachment Abutment	P5-3.75,4
392	Ball Attachment Abutment	P5-3.75,5
393	Ball Attachment Abutment	P5-3.75,6
394	Ball Attachment for I6B Implant	P5-I6B-1
395	Ball Attachment for I6B Implant	P5-I6B-2
396	Ball Attachment for I6B Implant	P5-I6B-3
397	Ball Attachment for I6B Implant -a+b+c	P5-I6B-1set
398	Ball Attachment for I6B Implant -a+b+c	P5-I6B-2set
399	Ball Attachment for I6B Implant - a+b+c	P5-I6B-3set
400	Angular Ball Attachment 20° - 1 mm - a+ b+c	P5-3.75,20-1set
401	Angular Ball Attachment 20° - 2 mm - a+ b+c	P5-3.75,20-2set
402	Angular Ball Attachment 20° - 3 mm - a+ b+c	P5-3.75,20-3set
403	Angular Ball Attachment 20° - 4 mm - a+ b+c	P5-3.75,20-4set
404	Angular Ball Attachment 20° - 5 mm - a+ b+c	P5-3.75,20-5set
405	Angular Ball Attachment 20° - 6 mm - a+ b+c	P5-3.75,20-6set
406	Angular Ball Attachment 20° - 1 mm	P5-3.75,20-1
407	Angular Ball Attachment 20° - 2 mm	P5-3.75,20-2
408	Angular Ball Attachment 20° - 3 mm	P5-3.75,20-3
409	Angular Ball Attachment 20° - 4 mm	P5-3.75,20-4
410	Angular Ball Attachment 20° - 5 mm	P5 -3.75,20-5
411	Angular Ball Attachment 20° - 6 mm	P5 -3.75,20-6
412	Ball Attachment - 1 mm	P5-4.5,1
413	Ball Attachment - 2 mm	P5-4.5,2
414	Ball Attachment - 4 mm	P5-4.5,4
415	Ball Attachment	P5-4.5,1set
416	Ball Attachment	P5-4.5,2set
417	Ball Attachment	P5-4.5,4set
418	Ball for Angular Adaptor	P5-P14,1
419	Ball for Angular Adaptor	P5-P14,2
420	Aesthetic Abutment with Hex	P6H-3.75,0.5
421	Aesthetic Abutment with Hex	P6H-3.75,1.5
422	Aesthetic Abutment with Hex	P6H-3.75,3
423	Aesthetic Abutment with Hex	P6H-3.75,4
424	Aesthetic Abutment with Hex	P6H-3.75,6

425	Aesthetic Abutment with Hex	P6H-4.5,1
426	Aesthetic Abutment with Hex	P6H-4.5,2
427	Aesthetic Abutment with Hex	P6H-4.5,3
428	Titanium Sleeve without Hex	P6HbT
429	Plastic Sleeve without Hex	P6H-b
430	Plastic Sleeve with Hex	P6H-bH
431	Titanium Sleeve with Hex	P6HbT-H
432	Screw	P6Ha
433	Screw for P6	P6a
434	Screw for P6- Long	P6a-L
435	Anti-rotation Aesthetic Abutment	P6-3.75,3
436	Anti-rotation Aesthetic Abutment	P6-3.75,4
437	Anti-rotation Aesthetic Abutment	P6-3.75,6
438	Anti-rotation Aesthetic Abutment	P6-4.5,3
439	Anti-rotation Aesthetic Abutment	P6-4.5,4
440	Anti-rotation Aesthetic Abutment	P6-4.5,6
441	Anti-rotation Aesthetic Abutment	P7-3.75,1
442	Anti-rotation Aesthetic Abutment	P7-3.75,2
443	Anti-rotation Aesthetic Abutment	P7-3.75,3
444	Plastic Sleeve without Hex	P7b
445	Plastic Sleeve with hex	P7b-H
446	Screw	P7a-1
447	Screw	P7a-2
448	Screw	P7a-3
449	Plastic abutment for casting	P9-3.75,11
450	Plastic abutment for casting - Perfect Fit	P9-3.75,11-P
451	Plastic abutment for casting	P9-4.5,11
452	Plastic abutment for casting - Perfect Fit	P9-4.5,11-P
453	GOLD Composed NOT Hex Abutment	P9G-3.75
454	GOLD Composed NOT Hex Abutment	P9G-4.5
455	Composed Hex Abutment	P9H-3.75,11
456	GOLD Composed NOT Hex Abutment	P9G-3.75,11
457	GOLD Composed Hex Abutment	P9HG-3.75,11
458	GOLD Composed Hex Abutment	P9HG-4.5
459	Flat Connection Abutment (base + sleeve + screw)	P12-3.75
460	Temporary Flat Connection Abutment	P12-3.75-T
461	Temporary Flat Connection Abutment	P12-3.75-T/L
462	Adaptor for flat connection abutment	P12C-3.75
463	Adaptor for flat connection abutment	P12C-3.75-5
464	Screw for P12	P12a
465	Screw for P12C	P12Ca
466	Angular Adaptor	P14-3.75,17-1
467	Angular Adaptor	P14-3.75,17-3
468	Angular Adaptor	P14-3.75,30-1
469	Angular Adaptor	P14-3.75,30-3
470	Plastic Sleeve for Angular Adaptor	P14b
471	Titanium Sleeve for Angular Adaptor	P14-b/T
472	Gold Sleeve for Angular Adaptor	P14-b/G
473	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3,1
474	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3,2
475	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3,3
476	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3,4
477	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3,5

478	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3,6
479	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3,7
480	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3.75,1
481	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3.75,2
482	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3.75,3
483	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3.75,4
484	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3.75,5
485	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3.75,6
486	Temporary Abutment for Immediate Loading (abutment + plastic cap)	P15-3.75,7
487	Plastic Cap for Temporary Abutment for Immediate Loading	P15a-1
488	Plastic Cap for Temporary Abutment for Immediate Loading	P15a-2
489	Plastic Cap for Temporary Abutment for Immediate Loading	P15a-3
490	Straight Adaptor	P16-3.75,1
491	Straight Adaptor	P16-3.75,2
492	Straight Adaptor	P16-3.75,3
493	Straight Adaptor	P16-3.75,4
494	Straight Adaptor	P16-3.75,5
495	Straight Adaptor	P16-3,1
496	Straight Adaptor	P16-3,2
497	Straight Adaptor	P16-3,3
498	Straight Adaptor	P16-3,4
499	Straight Adaptor	P16-3,5
500	Angular Adaptor	P24-3.75,17-1
501	Angular Adaptor	P24-3.75,17-3
502	Angular Adaptor	P24-3.75,30-1
503	Angular Adaptor	P24-3.75,30-3
504	Angular Adaptor	P24-3.75,45-1
505	Angular Adaptor	P24-3.75,45-3
506	Locator® Diameter 3.75, Length 0	P25-3.75,0
507	Locator® Diameter 3.75, Length 1	P25-3.75,1
508	Locator® Diameter 3.75, Length 2	P25-3.75,2
509	Locator® Diameter 3.75, Length 3	P25-3.75,3
510	Locator® Diameter 3.75, Length 4	P25-3.75,4
511	Locator® Diameter 3.75, Length 5	P25-3.75,5
512	Locator® Diameter 3.75, Length 6	P25-3.75,6
513	Locator® Diameter 3.75, Length 7	P25-3.75,7
514	Locator® Diameter 3.75, Length 8	P25-3.75,8
515	Locator® Diameter 3.75, Length 9	P25-3.75,9
516	Metal cover	P25a
517	Locator® Male Processing Package Assembly	P25-a,b/10
518	Locator® Extended Range Male Processing Package	P25-a,b/20
519	AB LOC Width 3.75, length 0 mm	P35-3.75,0
520	AB LOC Width 3.75, length 1 mm	P35-3.75,1
521	AB LOC Width 3.75, length 2 mm	P35-3.75,2
522	AB LOC Width 3.75, length 3 mm	P35-3.75,3
523	AB LOC Width 3.75, length 4 mm	P35-3.75,4
524	AB LOC Width 3.75, length 5 mm	P35-3.75,5
525	AB LOC Width 3.75, length 6 mm	P35-3.75,6
526	AB LOC Width 3.75, length 7 mm	P35-3.75,7
527	AB LOC Width 3.75, length 8 mm	P35-3.75,8
528	Low Connector Diameter 3, Length 1	P55-3,1
529	Low Connector Diameter 3, Length 2	P55-3,2

530	Low Connector Diameter 3, Length 3	P55-3,3
531	Low Connector Diameter 3, Length 4	P55-3,4
532	Low Connector Diameter 3.75, Length 1	P55-3.75,1
533	Low Connector Diameter 3.75, Length 2	P55-3.75,2
534	Low Connector Diameter 3.75, Length 3	P55-3.75,3
535	Low Connector Diameter 3.75, Length 4	P55-3.75,4
536	Low Connector Diameter 3.75, Length 5	P55-3.75,5
537	Low Connector Diameter 3.75, Length 6	P55-3.75,6
538	Low Connector Diameter 3.75, Length 7	P55-3.75,7
539	Low Connector Diameter 3.75, Length 8	P55-3.75,8
540	Metal Cap	P55a
541	Silicon Cap	P55b-1
542	Silicon Cap	P55b-2
543	Silicon Cap	P55b-3
544	Analog	D1-3
545	Analog	D1-3.75
546	Analog	D1-4.5
547	Analog	D1-5
548	Analog	D1-6
549	Implant Analog	D1-I6
550	Implant Analog	D1-I6B
551	Implant Analog	D1-P5-I6B
552	Implant Analog for P6	D1-P6
553	Implant Analog for P7	D1-P7
554	Analog for Aesthetic Abutment	D1-P6H
555	Analog for P14	D1-P14
556	Impression Transfer	D1K-3.75
557	Impression Transfer	D2-3,9
558	Impression Transfer	D2-3,15
559	Impression Transfer for open tray	D2O-3.75,9
560	Plastic Transfer	D2K-3.75
561	Impression Transfer for open tray	D2O-4.5,9
562	Impression Transfer for open tray	D2O-4.5,15
563	Impression Transfer for closed tray	D2-4.5,9
564	Impression Transfer for closed tray	D2-4.5,15
565	Narrow Impression Transfer for open tray	D2NO-3.75,9
566	Narrow Impression Transfer for closed tray	D2N-3.75,9
567	Impression Transfer for open tray	D2O-3.75,15
568	Narrow Impression Transfer for open tray	D2NO-3.75,15
569	Narrow Impression Transfer for closed tray	D2N-3.75,15
570	Narrow Transfer for P12	D2N-P12-3.75,9
571	Narrow Transfer for P12	D2N-P12-3.75,15
572	Impression Transfer for I6 Implant	D2-I6
573	Impression Transfer for I6B Implant	D2-I6B
574	Transfer for I7 diameter 3.75	D2-I7-3.75
575	Transfer for I7 diameter 4.2	D2-I7-4.2
576	Transfer for I7 diameter 5	D2-I7-5
577	Transfer for Aesthetic Abutment	D2-P6H
578	Transfer for Flat Connection Abutment	D2-P12-3.75,9
579	Transfer for Flat Connection Abutment	D2-P12-3.75,15
580	Transfer for Flat Connection Abutment	D2-P12-4.5,9
581	Transfer for Angular Adaptor	D2-P14

582	Impression Transfer for P6	D2-P6
583	Transfer for P7	D2-P7
584	Impression Transfer with Short screw	D2-3.75,9
585	Impression Transfer with Short screw	D2-3.75,15
586	Short fixing screw	D2a
587	Long fixing screw	D2a1
588	Clip Transfer	D3-3.75,9
589	Clip Transfer	D3-3.75,15
590	Narrow Clip Transfer	D3N-3.75,9
591	Narrow Clip Transfer	D3N-3.75,15
592	Clip Transfer	D3-4.5,9
593	Clip Transfer	D3-4.5,15
594	Prosthetic Kit	PK-3.75,1
595	Healing Cap	PK - Po
596	Plastic Conical Sleeve, (Red)	PK-P2
597	Plastic Conical Sleeve, (White)	PK-P2H
598	Anatomic Anti -rotation Abutment	PK-P3-1
599	Implant Analog	PK-D1
600	Plastic Transfer	PK-D2
601	Prosthetic Kit	PK-3.75,2
602	Anatomic Anti -rotation Abutment	PK-P3-2
603	Prosthetic Kit	PK-3.75,3
604	Anatomic Anti -rotation Abutment	PK-P3-3
605	Prosthetic Kit	PK-3.75,4
606	Anatomic Anti -rotation Abutment	PK-P3-4
607	Prosthetic Kit	PK-3,1
608	Prosthetic Kit	PK-3,2
609	Prosthetic Kit	PK-3,3
610	Prosthetic Kit	PK-3,4
611	Ratchet Hex Driver for Abutment	T1-1.2,9
612	Ratchet Hex Driver for Abutment	T1-1.2,15
613	Ratchet Hex Driver for P6	T1-1.8,9
614	Ratchet Hex Driver for P6	T1-1.8,15
615	Hand Hex Driver for Abutment	T2-1.2,9
616	Hand Hex Driver for Abutment	T2-1.2,15
617	Hand Hex Driver for P6	T2-1.8,9
618	Hand Hex Driver for P6	T2-1.8,15
619	T2-Drill Stopper Kit	T2-DDL
620	Hand Hex Driver for P7	T2-P7
621	Manual Adaptor for Driver	T3a
622	Ratchet Hex Driver for Implant	T3-2,9
623	Ratchet Hex Driver for Implant	T3-2,18
624	Ratchet Hex Driver for Implant	T3-2.4,9
625	Ratchet Hex Driver for Implant	T3-2.4,15
626	Ratchet Hex Driver for Implant	T3-2.4,18
627	Ratchet Hex Driver for I6B	T3-I6B-S
628	Ratchet Hex Driver for I6B	T3-I6B-L
629	Ratchet Hex Driver for I6B	T3-I6B-18
630	Ratchet Hex Driver for Implant	T3-I6-L
631	Retrieving Screw	T4
632	Contra Angle Driver for implant	T5-2.4,20
633	Contra Angle Driver for implant	T5F-2.4,20

634	Contra Angle Driver for implant	T5-2.4,25
635	Contra Angle Driver for implant	T5F-2.4,25
636	Contra Angle Driver for abutment	T5-1.2,21
637	Contra Angle Driver for abutment	T5-1.2,26
638	Flexible Driver for Hand piece	T5-1.2F
639	Contra Angle Driver for implant	T5-2,20
640	Contra Angle Driver for implant	T5F-2,20
641	Contra Angle Driver for implant	T5-2,25
642	Contra Angle Driver for implant	T5F-2,25
643	Ratchet Wrench	T8
644	Ratchet-Torque Combination	T8c-10-40
645	Depth Gauge	T9
646	Handle	T10
647	Mallet	T11
648	Abutment Gripper	T12
649	Technician's Handle	T13
650	Tap-3.25	T14-3.25
651	Tap-3.75	T14-3.75
652	Tap-4.2	T14-4.2
653	Tap-5	T14-5
654	Tissue Punch Driver	T17
655	Adaptor ratchet for Locator®	T1-T25
656	LOCATOR® Loc- Square Drive Torque Wrench Driver -15 mm	T25-c/15
657	LOCATOR® Loc- Square Drive Torque Wrench Driver -21 mm	T25-c/21
658	LOCATOR® Core Tool	T25-b
659	Parallel Pin	TP-17
660	Parallel Pin	TP-23
661	Ø1.9 Marker Drill Bit	TMD -1.9
662	Stopper Sleeve for drill	TS-2,5
663	Stopper Sleeve for drill	TS-2.5,5
664	Stopper Sleeve for drill	TS-2.8,5
665	Stopper Sleeve for drill	TS-3.2,5
666	Stopper Sleeve for drill	TS-3.65,5
667	Stopper Sleeve for drill	TS-4,5
668	Stopper Sleeve for drill	TS-4.5,5
669	Stopper Sleeve for drill	TS-5.5,5
670	Ø 2 Pilot Drill Bit	TPD -2
671	Ø 2.5 Drill Bit	TD-2.5
672	Ø 2.8 Drill Bit	TD-2.8
673	Ø 3.2 Drill Bit	TD-3.2
674	Ø 3.65 Drill Bit	TD-3.65
675	Ø 3.8 Drill Bit	TD-3.8
676	Ø 4.0 Drill Bit	TD-4.0
677	Ø 4.2 Drill Bit	TD-4.2
678	Ø 4.5 Drill Bit	TD-4.5
679	Ø 5Drill Bit	TD-5
680	Ø 5.2 Drill Bit	TD-5.2
681	Ø 5.5 Drill Bit	TD-5.5
682	Drill extension Bit with irrigation	TDE-I
683	Ø1.9 Marker Drill Bit with irrigation	TMDI -1.9
684	Ø 2 Pilot Drill Bit with irrigation	TPDI -2.0
685	Ø 2 Pilot Conical Drill Bit	TPDC-2.0

686	ø 2 Pilot Conical Drill Bit with irrigation	TPDCI-2.0
687	ø 2.5 Drill Bit with irrigation	TDI-2.5
688	ø 2.8 Drill Bit with irrigation	TDI-2.8
689	ø 3.2 Drill Bit with irrigation	TDI-3.2
690	ø 3.65 Drill Bit with irrigation	TDI-3.65
691	ø 3.8 Drill Bit with irrigation	TDI-3.8
692	ø 4.0 Drill Bit with irrigation	TDI-4.0
693	ø 4.2 Drill Bit with irrigation	TDI-4.2
694	ø 4.5 Drill Bit with irrigation	TDI-4.5
695	ø 5 Drill Bit with irrigation	TDI-5
696	ø 5.2 Drill Bit with irrigation	TDI-5.2
697	ø 5.5 Drill Bit with irrigation	TDI-5.5
698	ø 2.2 Conical Drill Bit	TDC-2.2
699	ø 2.7 Conical Drill Bit	TDC-2.7
700	ø 3.3 Conical Drill Bit	TDC-3.3
701	ø 3.7 Conical Drill Bit	TDC-3.7
702	ø 4.0 Conical Drill Bit	TDC-4.0
703	ø 4.2 Conical Drill Bit	TDC-4.2
704	ø 4.5 Conical Drill Bit	TDC-4.5
705	ø 5.5 Conical Drill Bit	TDC-5.5
706	Counter Sink	TDCS-3.7-4.2
707	Counter Sink	TDCS-5-6
708	Conical Stopper Drill	TDCSt-2.2,10
709	Conical Stopper Drill	TDCSt-2.2,11.5
710	Conical Stopper Drill	TDCSt-2.2,13
711	Conical Stopper Drill	TDCSt-2.2,6
712	Conical Stopper Drill	TDCSt-2.2,8
713	Conical Stopper Drill	TDCSt-2.7,10
714	Conical Stopper Drill	TDCSt-2.7,11.5
715	Conical Stopper Drill	TDCSt-2.7,13
716	Conical Stopper Drill	TDCSt-2.7,6
717	Conical Stopper Drill	TDCSt-2.7,8
718	Conical Stopper Drill	TDCSt-3.3,10
719	Conical Stopper Drill	TDCSt-3.3,11.5
720	Conical Stopper Drill	TDCSt-3.3,13
721	Conical Stopper Drill	TDCSt-3.3,6
722	Conical Stopper Drill	TDCSt-3.3,8
723	Conical Stopper Drill	TDCSt-3.7,10
724	Conical Stopper Drill	TDCSt-3.7,11.5
725	Conical Stopper Drill	TDCSt-3.7,13
726	Conical Stopper Drill	TDCSt-3.7,6
727	Conical Stopper Drill	TDCSt-3.7,8
728	Conical Stopper Drill	TDCSt-4.0,10
729	Conical Stopper Drill	TDCSt-4.0,11.5
730	Conical Stopper Drill	TDCSt-4.0,13
731	Conical Stopper Drill	TDCSt-4.0,6
732	Conical Stopper Drill	TDCSt-4.0,8
733	Conical Stopper Drill	TDCSt-4.5,10
734	Conical Stopper Drill	TDCSt-4.5,11.5
735	Conical Stopper Drill	TDCSt-4.5,13
736	Conical Stopper Drill	TDCSt-4.5,6
737	Conical Stopper Drill	TDCSt-4.5,8

738	Conical Stopper Drill	TDCSt-5.5,10
739	Conical Stopper Drill	TDCSt-5.5,11.5
740	Conical Stopper Drill	TDCSt-5.5,13
741	Conical Stopper Drill	TDCSt-5.5,6
742	Conical Stopper Drill	TDCSt-5.5,8
743	Conical Stopper Drill with irrigation	TDCStI-2.2,10
744	Conical Stopper Drill with irrigation	TDCStI-2.2,11.5
745	Conical Stopper Drill with irrigation	TDCStI-2.2,13
746	Conical Stopper Drill with irrigation	TDCStI-2.2,6
747	Conical Stopper Drill with irrigation	TDCStI-2.2,8
748	Conical Stopper Drill with irrigation	TDCStI-2.7,10
749	Conical Stopper Drill with irrigation	TDCStI-2.7,11.5
750	Conical Stopper Drill with irrigation	TDCStI-2.7,13
751	Conical Stopper Drill with irrigation	TDCStI-2.7,6
752	Conical Stopper Drill with irrigation	TDCStI-2.7,8
753	Conical Stopper Drill with irrigation	TDCStI-3.3,10
754	Conical Stopper Drill with irrigation	TDCStI-3.3,11.5
755	Conical Stopper Drill with irrigation	TDCStI-3.3,13
756	Conical Stopper Drill with irrigation	TDCStI-3.3,6
757	Conical Stopper Drill with irrigation	TDCStI-3.3,8
758	Conical Stopper Drill with irrigation	TDCStI-3.7,10
759	Conical Stopper Drill with irrigation	TDCStI-3.7,11.5
760	Conical Stopper Drill with irrigation	TDCStI-3.7,13
761	Conical Stopper Drill with irrigation	TDCStI-3.7,6
762	Conical Stopper Drill with irrigation	TDCStI-3.7,8
763	Conical Stopper Drill with irrigation	TDCStI-4.0,10
764	Conical Stopper Drill with irrigation	TDCStI-4.0,11.5
765	Conical Stopper Drill with irrigation	TDCStI-4.0,13
766	Conical Stopper Drill with irrigation	TDCStI-4.0,6
767	Conical Stopper Drill with irrigation	TDCStI-4.0,8
768	Conical Stopper Drill with irrigation	TDCStI-4.5,10
769	Conical Stopper Drill with irrigation	TDCStI-4.5,11.5
770	Conical Stopper Drill with irrigation	TDCStI-4.5,13
771	Conical Stopper Drill with irrigation	TDCStI-4.5,6
772	Conical Stopper Drill with irrigation	TDCStI-4.5,8
773	Conical Stopper Drill with irrigation	TDCStI-5.5,10
774	Conical Stopper Drill with irrigation	TDCStI-5.5,11.5
775	Conical Stopper Drill with irrigation	TDCStI-5.5,13
776	Conical Stopper Drill with irrigation	TDCStI-5.5,6
777	Conical Stopper Drill with irrigation	TDCStI-5.5,8
778	Drill Bit with irrigation	TDI-I15-4.5
779	Drill Bit with irrigation	TDI-I15-5.5
780	Drill Bit with irrigation	TDI-I15-6.5
781	Drill Bit with irrigation	TDI-I15-7.5
782	Drill Bit with irrigation	TDI-I15-8.5
783	∅ 2.2mm Conical Drill Bit with irrigation	TDCI-2.2
784	∅ 2.4mm Conical Drill Bit with irrigation	TDCI-2.4
785	∅ 2.7mm Conical Drill Bit with irrigation	TDCI-2.7
786	∅ 3.3mm Conical Drill Bit with irrigation	TDCI-3.3
787	∅ 3.65mm Conical Drill Bit with irrigation	TDCI-3.65
788	∅ 3.7mm Conical Drill Bit with irrigation	TDCI-3.7
789	∅ 4.0mm Conical Drill Bit with irrigation	TDCI-4.0

790	ø 2.2mm Conical Drill Bit with irrigation	TDCI -4.2
791	ø 4.5mm Conical Drill Bit with irrigation	TDCI -4.5
792	ø 5.2mm Conical Drill Bit with irrigation	TDCI -5.2
793	ø 5.5mm Conical Drill Bit with irrigation	TDCI -5.5
794	ø6.5mm,Length10mm Conical Drill Bit with irrigation	TDCI -6.5,10
795	ø6.5mm,Length11.5mm Conical Drill Bit with irrigation	TDCI -6.5,11.5
796	ø 6.5mm Length,8 mm Conical Drill Bit with irrigation	TDCI -6.5,8
797	ø 7.5mm, Length 10mm Conical Drill Bit with irrigation	TDCI -7.5,10
798	ø 7.5mm, Length 11.5mm Conical Drill Bit with irrigation	TDCI -7.5,11.5
799	ø7mm. Length 5,8mm Conical Drill Bit with irrigation	TDCI -7.5,8
800	Stopper Drill Bit Diameter 2mm Length 10mm	TDSt-2,10
801	Stopper Drill Bit Diameter 2mm Length 11.5mm	TDSt-2,11.5
802	Stopper Drill Bit Diameter 2mm Length 13mm	TDSt-2,13
803	Guid Drill Bit Diameter 2mm Length 17mm	TDG-2,17
804	Guid Drill Bit Diameter 2mm Length 21mm	TDG-2,21
805	Guid Drill Bit Diameter 2mm Length 25mm	TDG-2,25
806	Guid Drill Bit Diameter 2mm Length 30mm	TDG-2,30
807	Stopper Drill Bit Diameter 2mm Length 6mm	TDSt-2,6
808	Stopper Drill Bit Diameter 2mm Length 8mm	TDSt-2,8
809	Stopper Drill Bit Diameter 2.5mm Length 6mm	TDSt-2.5,6
810	Stopper Drill Bit Diameter 2.5mm Length 8mm	TDSt-2.5,8
811	Stopper Drill Bit Diameter 2.5mm Length 10mm	TDSt-2.5,10
812	Stopper Drill Bit Diameter 2.5mm Length 11.5mm	TDSt-2.5,11.5
813	Stopper Drill Bit Diameter 2.5mm Length 13mm	TDSt-2.5,13
814	Guid Drill Bit Diameter 2.5mm Length 17mm	TDG-2.5,17
815	Guid Drill Bit Diameter 2.5mm Length 21mm	TDG-2.5,21
816	Guid Drill Bit Diameter 2.5mm Length 25mm	TDG-2.5,25
817	Guid Drill Bit Diameter 2.5mm Length 30mm	TDG-2.5,30
818	Stopper Drill Bit Diameter 2.8mm Length 6mm	TDSt-2.8,6
819	Stopper Drill Bit Diameter 2.8mm Length 8mm	TDSt-2.8,8
820	Stopper Drill Bit Diameter 2.8mm Length 10mm	TDSt-2.8,10
821	Stopper Drill Bit Diameter 2.8mm Length 11.5mm	TDSt-2.8,11.5
822	Stopper Drill Bit Diameter 2.8mm Length 13mm	TDSt-2.8,13
823	Guid Drill Bit Diameter 2.8mm Length 17mm	TDG-2.8,17
824	Guid Drill Bit Diameter 2.8mm Length 21mm	TDG-2.8,21
825	Guid Drill Bit Diameter 2.8mm Length 25mm	TDG-2.8,25
826	Guid Drill Bit Diameter 2.8mm Length 30mm	TDG-2.8,30
827	Stopper Drill Bit Diameter 3.2mm Length 6mm	TDSt-3.2,6
828	Stopper Drill Bit Diameter 3.2mm Length 8mm	TDSt-3.2,8
829	Stopper Drill Bit Diameter 3.2mm Length 10mm	TDSt-3.2,10
830	Stopper Drill Bit Diameter 3.2mm Length 11.5mm	TDSt-3.2,11.5
831	Stopper Drill Bit Diameter 3.2mm Length 13mm	TDSt-3.2,13
832	Guid Drill Bit Diameter 3.2mm Length 17mm	TDG-3.2,17
833	Guid Drill Bit Diameter 3.2mm Length 21mm	TDG-3.2,21
834	Guid Drill Bit Diameter 3.2mm Length 25mm	TDG-3.2,25
835	Guid Drill Bit Diameter 3.2mm Length 30mm	TDG-3.2,30
836	Stopper Drill Bit Diameter 3.65mm Length 6mm	TDSt-3.65,6
837	Stopper Drill Bit Diameter 3.65mm Length 8mm	TDSt-3.65,8
838	Stopper Drill Bit Diameter 3.65mm Length 10mm	TDSt-3.65,10
839	Stopper Drill Bit Diameter 3.65mm Length 13mm	TDSt-3.65,11.5
840	Stopper Drill Bit Diameter 3.65mm Length 25mm	TDSt-3.65,13
841	Guid Drill Bit Diameter 3.65mm Length 17mm	TDG-3.65,17

842	Guid Drill Bit Diameter 3.65mm Length 21mm	TDG-3.65,21
843	Guid Drill Bit Diameter 3.65mm Length 25mm	TDG-3.65,25
844	Guid Drill Bit Diameter 3.65mm Length 30mm	TDG-3.65,30
845	Stopper Drill Bit Diameter 4mm Length 6mm	TDSt-4,6
846	Stopper Drill Bit Diameter 4mm Length 8mm	TDSt-4,8
847	Stopper Drill Bit Diameter 4mm Length 10mm	TDSt-4,10
848	Stopper Drill Bit Diameter 4mm Length 11.5mm	TDSt-4,11.5
849	Stopper Drill Bit Diameter 4mm Length 13mm	TDSt-4,13
850	Guid Drill Bit Diameter 4mm Length 17mm	TDG-4,17
851	Guid Drill Bit Diameter 4mm Length 21mm	TDG-4,21
852	Guid Drill Bit Diameter 4mm Length 25mm	TDG-4,25
853	Stopper Drill Bit Diameter 4.2mm Length 6mm	TDSt-4.2,6
854	Stopper Drill Bit Diameter 4.2mm Length 8mm	TDSt-4.2,8
855	Stopper Drill Bit Diameter 4.2mm Length 10mm	TDSt-4.2,10
856	Stopper Drill Bit Diameter 4.2mm Length 11.5mm	TDSt-4.2,11.5
857	Stopper Drill Bit Diameter 4.2mm Length 13mm	TDSt-4.2,13
858	Stopper Drill Bit Diameter 4.5mm Length 6mm	TDSt-4.5,6
859	Stopper Drill Bit Diameter 4.5mm Length 8mm	TDSt-4.5,8
860	Stopper Drill Bit Diameter 4.5mm Length 10mm	TDSt-4.5,10
861	Stopper Drill Bit Diameter 4.5mm Length 11.5mm	TDSt-4.5,11.5
862	Stopper Drill Bit Diameter 4.5mm Length 13mm	TDSt-4.5,13
863	Guid Drill Bit Diameter 4.5mm Length 17mm	TDG-4.5,17
864	Guid Drill Bit Diameter 4.5mm Length 21mm	TDG-4.5,21
865	Guid Drill Bit Diameter 4.5mm Length 25mm	TDG-4.5,25
866	Guid Drill Bit Diameter 5.5mm Length 17mm	TDG-5.5,17
867	Guid Drill Bit Diameter 5.5mm Length 21mm	TDG-5.5,21
868	Guid Drill Bit Diameter 5.5mm Length 25mm	TDG-5.5,25
869	Punch for cutting gums	TD-T17
870	Punch for cutting gums	TD-T17-4.35
871	Trephine	TDT-3.0
872	Trephine	TDT-4.0
873	Trephine	TDT-5.0
874	Drilling Extension Bit	TD-E
875	Short Drill	TPDS-2.0
876	Short Drill	TDS-2.5
877	Short Drill	TDS-2.8
878	Short Drill	TDS-3.2
879	Short Drill	TDS-3.65
880	Short Drill	TDS-4.2
881	Short Drill	TDS-5.2
882	Short Drill with irrigation	TPDSI-2.0
883	Stopper Drill Bit with irrigation	TDStI-2,10
884	Stopper Drill Bit with irrigation	TDStI-2,11.5
885	Stopper Drill Bit with irrigation	TDStI-2,13
886	Stopper Drill Bit with irrigation	TDStI-2,6
887	Stopper Drill Bit with irrigation	TDStI-2,8
888	Short Drill with irrigation	TDSI-2.0
889	Short Drill with irrigation	TDSI-2.5
890	Stopper Drill Bit with irrigation	TDStI-2.5,10
891	Stopper Drill Bit with irrigation	TDStI-2.5,11.5
892	Stopper Drill Bit with irrigation	TDStI-2.5,13
893	Stopper Drill Bit with irrigation	TDStI-2.5,6
894	Stopper Drill Bit with irrigation	TDStI-2.5,8

895	Short Drill with irrigation	TDSI-2.8
896	Stopper Drill Bit with irrigation	TDStI-2.8,10
897	Stopper Drill Bit with irrigation	TDStI-2.8,11.5
898	Stopper Drill Bit with irrigation	TDStI-2.8,13
899	Stopper Drill Bit with irrigation	TDStI-2.8,6
900	Stopper Drill Bit with irrigation	TDStI-2.8,8
901	Short Drill with irrigation	TDSI-3.2
902	Stopper Drill Bit with irrigation	TDStI-3.2,6
903	Stopper Drill Bit with irrigation	TDStI-3.2,8
904	Stopper Drill Bit with irrigation	TDStI-3.2,10
905	Stopper Drill Bit with irrigation	TDStI-3.2,11.5
906	Stopper Drill Bit with irrigation	TDStI-3.2,13
907	Short Drill with cooling	TDSI-3.65
908	Stopper Drill Bit with irrigation	TDStI-3.65,6
909	Stopper Drill Bit with irrigation	TDStI-3.65,8
910	Stopper Drill Bit with irrigation	TDStI-3.65,10
911	Stopper Drill Bit with irrigation	TDStI-3.65,11.5
912	Stopper Drill Bit with irrigation	TDStI-3.65,13
913	Short Drill with irrigation	TDSI-4.0
914	Stopper Drill Bit with irrigation	TDStI-4,6
915	Stopper Drill Bit with irrigation	TDStI-4,8
916	Stopper Drill Bit with irrigation	TDStI-4,10
917	Stopper Drill Bit with irrigation	TDStI-4,11.5
918	Stopper Drill Bit with irrigation	TDStI-4,13
919	Short Drill with irrigation	TDSI-4.2
920	Stopper Drill Bit with irrigation	TDStI -4.2,6
921	Stopper Drill Bit with irrigation	TDStI -4.2,8
922	Stopper Drill Bit with irrigation	TDStI -4.2,10
923	Stopper Drill Bit with irrigation	TDStI -4.2,11.5
924	Stopper Drill Bit with irrigation	TDStI -4.2,13
925	Short Drill with irrigation	TDSI-4.5
926	Stopper Drill Bit with irrigation	TDStI -4.5,6
927	Stopper Drill Bit with irrigation	TDStI -4.5,8
928	Stopper Drill Bit with irrigation	TDStI -4.5,10
929	Stopper Drill Bit with irrigation	TDStI -4.5,11.5
930	Stopper Drill Bit with irrigation	TDStI -4.5,13
931	Short Drill with irrigation	TDSI-5.2
932	Short Drill with irrigation	TDSI-5.5
933	Long Drill	TDL-2.5
934	Trephine with irrigation	TDTI-3.0
935	Trephine with irrigation	TDTI-4.0
936	Trephine with irrigation	TDTI-5.0
937	Long Drill	TDL-2.8
938	Long Drill	TDL-3.2
939	Long Drill	TDL-3.65
940	Long Drill	TDL-4.2
941	Long Drill with irrigation	TDLI-2.5
942	Long Drill with irrigation	TDLI-2.8
943	Long Drill with irrigation	TDLI-3.2
944	Long Drill with irrigation	TDLI-3.65
945	Long Drill with irrigation	TDLI-4.2
946	Handle Diameter 2.0mm	TH-2.0

947	Handle Diameter 2.5mm	TH-2.5
948	Handle Diameter 2.8mm	TH-2.8
949	Handle Diameter 3.2mm	TH-3.2
950	Handle Diameter 3.65mm	TH-3.65
951	Pin Fixation	TP-2/31
952	Zirconium Kit	AB Shine
953	Pointed head burs for inter proximal preparation	PR15
954	Versatile bur for all abutments aspects	K2
955	Versatile bur for all abutments aspects	PR2
956	Rounded head bur for preparation of the abutments top corners	A3
957	Rounded head bur for preparation of the abutments top corners	A4
958	Bur for abutment preparation above the anatomic level	T2
959	Bur for abutment preparation above the anatomic level	E5
960	Bur for abutment preparation above the anatomic level	E7
961	Sinus Lift Kit	TKSinus
962	Light Blue Sinus Lift Instrument	Tsinus-0
963	Green Sinus Lift Instrument	Tsinus-1
964	Red Sinus Lift Instrument	Tsinus-2
965	Black Sinus Lift Instrument	Tsinus-3
966	Silver Sinus Lift Instrument	Tsinus-4
967	Pink Sinus Lift Instrument	Tsinus-5
968	Purple Sinus Lift Instrument	Tsinus-6
969	Gold Sinus Lift Instrument	Tsinus-7
970	Straight Osteotom	TKOsteo
971	Straight Osteotom	T-ost-2-2.8
972	Straight Osteotom	T-ost-2.8-3.2
973	Straight Osteotom	T-ost-3.2-3.6
974	Straight Osteotom	T-ost-3.6-4.3
975	Straight Osteotom	T-ost-4.3-4.8
976	Osteotoms Kit	TKAOsteo
977	Angular Osteotom	T-ost/A-2-2.8
978	Angular Osteotom	T-ost/A-2.8-3.2
979	Angular Osteotom	T-ost/A-3.2-3.6
980	Angular Osteotom	T-ost/A-3.6-4.3
981	Angular Osteotom	T-ost/A-4.3-4.8
982	Resin Cement for Implant Retained Crowns	ABCem
983	TKD - Guided Service Drill Kit	TKD Guided

[<< Previous](#)

Legend:

^ Suspended Medical Device.

Last updated on 01 July 2014 © 2014. Health Sciences Authority. All Rights Reserved.

Best viewed using Internet Explorer 8.0 and above Privacy Statement / Terms of Use / HSA Data Protection Policy / Rate Our Website