




Superior Implant Technology



I22

THE ULTIMATE
IMPLANT FOR
IMMEDIATE LOADING

I22

THE ULTIMATE IMPLANT FOR IMMEDIATE LOADING



Features & Benefits:

- Conical implant with gentle threads
- Triple leads threads with distance of 0.9 between the threads
- Double thread at the apex, to provide better initial stability.
- Minimum pressure for type 1 bone.
- Recommend for bone type 1-2 lower jaw
- Triple lead to achieve intimate bone contact and shorter insertion time

Sizes: internal hex standard platform

Diameter (mm)	Platform	Length (mm)
3.75	Standard	8, 10, 11.5, 13, 16
4.2	Standard	8, 10, 11.5, 13, 16
5	Standard	8, 10, 11.5

www.ab-dent.com