




*Superior Implant Technology*

# I55


A NEW IMPLANT WITH  
THE OPTIMAL DESIGN FOR  
IMMEDIATE LOADING

---



# I55

## A NEW IMPLANT WITH THE OPTIMAL DESIGN FOR IMMEDIATE LOADING



### Features & Benefits:

- ▶ Spiral implant with combination of gentle and sharp threads for immediate loading:
  - Flat thread that enables the strongest initial stability
  - Sharp thread that enables the insertion of the implant easily and with minimal trauma to the bone
- ▶ A narrow apex for better primary stability
- ▶ Large surface area facilitating the osseointegration process (BIC)
- ▶ Available in narrow diameter
- ▶ Recommend for bone type 3-4 upper jaw

### Sizes: Standard & narrow platform, internal hexagon

Cat no.	Diameter	Platform	Length
I55	3	Narrow	10, 11.5, 13, 16
I55	3.3	Narrow	10, 11.5, 13, 16
I55	3.75	Standard	8, 10, 11.5, 13, 16
I55	4.2	Standard	8, 10, 11.5, 13, 16
I55	4.5	Standard	6, 8, 10, 11.5, 13, 16
I55	5	Standard	6, 8, 10, 11.5, 13

[www.ab-dent.com](http://www.ab-dent.com)